Sapphire Ceramic Kit

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Sapphire Ceramic Kit

Sapphire Ceramic Combi-Kit consists of :

Sapphire Ceramic, a light cure microhybrid composite with highest mechanical strength, suitable for all cavities. **Sapphire Ceramic** is based on BIS-GMA-resin and inorganic filler particles of 0.05-0.9 µm. Filler content: 81 % (wt), 63 % (vol.). **Sapphire Ceramic** meets the requirements of **DIN EN ISO 4049**.

- **Sapphire Conditioner**, a 37% phosphoric acid gel with an excellent thixotropy

- **Sapphire Bond**, a simple to use one bottle light cure single component priming and dentinoenamel adhesiv based on ethanol. It is designed for strong bonding of composites and metals to enamel and dentin and nonprecious and precious metals. **Sapphire Bond** can be also used for priming the root-canals before filling of the root canals or cementation of endodontic posts. Its strong adhesion to etched enamel or dentin proceeds on principles similar to that occurring with glassionomer cements. Good, long lasting adhesive strength and good biocompatibility are attained by methacrylate grafted polycarboxylic acids. Due to ist hydrophilicity, **Sapphire Bond** can be used on slightly moist dentin surfaces (wet-bonding-technique).

- **Sapphire Bond Activator**, a product for additional and/or extended applications of **Sapphire Bond** for all the special applications, where a dual or self cure priming and bonding system is strongly recommended.

- **Accessories**

Application:

1. Isolation

Rubber dam is the recommended method of isolation.

2. Cavity Preparation

Remove any debris and stain prior to preparation. Prepare the cavity with minimal tooth reduction. Margins should have a slight (**0.5 - 1.0 mm**) bevel placed in the enamel to increase the surface area for greater bond strength.

3. Pulp Protection

Cavity floor of deep excavations should be covered with a thin layer of hard setting calciumhydroxide material.

4. Enamel and Dentine Conditioning

Recommended is the total etch technique. Apply **Sapphire Bond** onto the enamel and dentine surfaces beginning with the enamel bevels. Leave the etching gel in place for **20 seconds**. On primary teeth and teeth high in fluoride, a **60 second** etch is recommended. Rinse for 20 seconds with water. Dry it in a water and oil-free airstream, but **do not** desiccate. A slightly wet dentine surface is important for the function of **Sapphire Bond**. The etched enamel bevel should have a chalky white appearance.

Etching precaution:

It is essential, that etched areas are not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water.

5. Application of Sapphire Bond

Apply **Sapphire Bond** generously with a brush onto the enamel and dentin surfaces for **30 seconds** with agitation. Dry cautiously with oil free air for 15 seconds. Do not desiccate the dentin.

Cure the **Sapphire Bond** coating by exposing its entire area to a dental halogen light unit for **20 seconds** before application of a second layer of **Sapphire Bond**.

Notes: Do not rinse off **Sapphire Bond**! If not used immediately, place dispensed **Sapphire Bond** in subdued light to prevent premature polymerisation by incident light. **Sapphire Bond** will **not** selfcure.

6. Application of a Second Layer of Sapphire Bond

Apply again **Sapphire Bond** generously with a brush onto the adhesive surfaces for **30 seconds** with agitation.

Thoroughly dry again the primed areas with oil free air for **15 seconds** and cure the **Sapphire Bond** coating by exposing its entire area to a dental halogen light unit for **20 seconds** before placement of a lightcure composite.

Instructions for Use

7. Colour Matching

Ascertain the tooth shade while teeth are still moist and select the appropriate **Sapphire Ceramic** shade.

Sapphire Ceramic is shaded according to VITA®-Shades.

8. Placement of Sapphire Ceramic

Apply **Sapphire Ceramic** in the selected shade. **Sapphire Ceramic** is easy to model. Place it by using a non-metallic instrument. Transparent matrix strips may be used.

9. Curing of Sapphire Ceramic

For an optimum result apply **Sapphire Ceramic** in **layers of max. 2 mm**. Light cure each increment separately.

By using a polymerization unit with a **light intensity** of **at least 1000mW/cm²**, cure each increment as follows:

Lighter shades (e.g. A1, A2, C2) **20 seconds**

Darker shades (e.g. A3.5, B3) **30 seconds**

10. Finishing

Contour with fine diamonds, stones or burs. Polish to high gloss with discs or rubber points. Interproximal finishing is accomplished by fine grit finishing strips.

11. Storage

Do not store above 25°C (78°F).

12. Additional Notes

- Do neither store the composite material in proximity of eugenol containing products, nor let the composite allow to come into contact with materials containing eugenol. Eugenol can impair the hardening of the composite and cause discoloration.
- Do not use any resin to adjust viscosity of composite restorative material.
- Contact of resin pastes with skin should be avoided, especially by anyone having known resin allergies.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.